



Emission Measurement System Designed to Meet Euro 7 and U.S. Tier 4 Standards

ECMA 17th International Conference ECT 2026

HORIBA, Ltd. Tomoshi Yoshimura

Oct 23, 2026



Agenda

1. Regulatory Trends
 - EURO 7 and U.S. Tier 4
2. Measurement Systems and Technical Challenges
3. Ultra-Low NO_x measurement
4. Low-Mass PM Measurement
5. Low-CO₂ with high moisture measurement
6. Summary

Euro 7, US Tier 4, and MY2027

Major Changes from the Previous Regulation

Regulation	Category	Date	Main Changes from the Previous Regulation	Requirement for Testing Equipment
Euro 7 	LDV 	Nov. 29 2026	<u>Updates</u> <ul style="list-style-type: none"> - PN23 to PN10 - PM and PN10 for all PI engines 	<ul style="list-style-type: none"> - PN10 Measurement
Euro 7 	HDV 	May 29, 2028	<u>New pollutants</u> <ul style="list-style-type: none"> - NMOG(Non-methane organic gas) - HCHO (TBD by 2027) - CH₄, N₂O <u>Updates</u> <ul style="list-style-type: none"> - NH₃ mass emission (Euro VI: ppm) - PN10 (Euro VI: PN23) 	<ul style="list-style-type: none"> - NMOG : Potential for the Impinger-Cartridge Method - HCHO : TBD - NH₃ Mass Measurement - CH₄, N₂O, PN10 Measurement
US Tier 4 	LDV 	MY2027 postponed to MY2029?	<u>Emission limit tightened:</u> <ul style="list-style-type: none"> - PM (3mg/mi → 0.5mg/mi) 	<ul style="list-style-type: none"> - Analytical Methods for Low-Mass PM measurement
MY2027 	HDV 	MY2027	<u>Emission limit tightened:</u> <ul style="list-style-type: none"> - NOx (200mg/hp-hr → 35mg/hp-hr) 	<ul style="list-style-type: none"> - Analytical Methods for Ultra-Low Concentrations NOx measurement

Euro 7 Chassis Test Cell Testing Configuration

LDV 



Euro 6e

Euro 7
Additional

Weather station



Driver's aid

*CH₄ required PI engine vehicle only**Gas Analysis**
CO/CO₂ / NO_x / THC / CH₄


Euro 7 HDD Engine Testing Configuration



Euro VI

Euro 7
Additional

*Whether HCHO measurement is required is TBD
*QCL-IR / FTIR analyzer with type approval must be used according to UNR49.




US Tier 4 NMOG + NOx, PM Emission Standards

LDV 

LDV


EPA (Environment Protection Agency) proposes to adjust the implementation timeline from model year 2027 to 2029 for the Tier4 criteria pollutant standards. (May 2026)

NMOG + NOx Fleet average emission limits



PM emission limits

2-years delayed, Phase-in MY2027 to MY2029 ?



Key items

Solution for Low-Mass PM measurements

US HD Regulation MY2027~

HD



Low NOx Standards will be implemented from MY2027.

HD

Current

Compression-Ignition Standards for Duty Cycle Testing

§ 86.007 11(a)

Duty Cycle	NOx mg/hp·hr	NMHC mg/hp·hr	PM mg/hp·hr	CO mg/hp·hr
FTP and SET	200	140	10	15.5

Spark-Ignition Standards for Duty Cycle Testing

§86.008 10(a)

Duty Cycle	NOx mg/hp·hr	NMHC mg/hp·hr	PM mg/hp·hr	CO mg/hp·hr
FTP and SET	200	140	10	14.4

NEW MY2027~

Compression-Ignition Standards for Duty Cycle Testing

§ 1036.104(1)

Duty Cycle	NOx mg/hp·hr	HC mg/hp·hr	PM mg/hp·hr	CO mg/hp·hr
SET and FTP	35	60	5	6.0
LLC	50	140	5	6.0

Spark-Ignition Standards for Duty Cycle Testing

§ 1036.104(2)

Duty Cycle	NOx mg/hp·hr	HC mg/hp·hr	PM mg/hp·hr	CO mg/hp·hr
SET and FTP	35	60	5	14.4
LLC	35	60	5	6.0

Key items

Solution for Ultra-Low NOx measurements

CO₂ Regulatory Trends


HDV and Urban Buses

EU

Latest regulation EU 2024/1610

Published May 14, 2024

CO₂ reduction requirement



- HDV : By 2030, a 45% reduction in CO₂ mass emissions from 2020 levels will be required.
- Urban buses : 100% reduction in emissions will be required from 2035 onward.

Note For Urban Buses, 100% reduction will be required from 2035.

NOx Emission Level

Reference Material : SwRI Research Result

Expected Low NOx Concentrations

Gas Condition : Conti-dilute
 Test Cycle : FTP Hot-start (HDDTC)
 Result : 0.010g/bhp-hr
Peak conc. : < 5ppm
Avg conc. (FWMC) : 0.19 ppm
Back ground : 0.055ppm

**This is the test case of SwRI.*

Some customers would get different results.

In many cases, minimum 10 ppm range is acceptable, but **2 ppm minimum range is recommended for low-NOx measurement where sub-1 ppm accuracy is necessary**



Expected concentration will be as follows:

- **Bag : less than 1ppm (Both cold-start and hot-start)**
- **Dilute-Conti : “meas” would be less than 1ppm (Both hot-start and cold-start)**

“meas” refers to the actual flow-weighted mean measured over any test interval (CFR 1065.205)

Concentration values for peak, FWMC, and background are very low compared to 2010 engines


Stage 3 Hot-Start at 0.010 g/hp-hr Dilute



Mass Result Impact Consideration

Hot FTP, Dilute Reference Material : SwRI data, reported in EMTC meeting


Raw vs Dilute Example Hot FTP – 0.019 g/hp-hr



- Stage I final matrix – Dilute and Raw in parallel
- Dilute at 0.0189 g/hp-hr
- Raw at 0.0203 g/hp-hr

SwRI Data


Stage 3 Hot-Start at 0.010 g/hp-hr Dilute



NOx mass rates keeps high level from 600-900 sec (continuous high rpm section) of hot-FTP phase.





The mass emission in this section account for much of the total mass emission, though concentration based emission keeps low level, **below 0.5 ppm**.

The analyzed concentration has a proportional impact on mass emission. More accurate analyzer at low concentration range will be helpful for more precise data.



Challenges in Ultra-Low-Concentration NOx Measurement

HD Tier 4, use case

	Cold-Start	Hot-Start
Dilution Gas	Highly accurate measurement of ultra-low-concentration NOx is required. Even slight drift or noise can have a major impact on calculated emissions mass.	
	 <p>NOx Concentration (ppm)</p> <p>Time (s)</p> <p>Average concentration: ~1 ppm</p> <p>Measurements in the ppb range are required; however, conventional standard analyzers are more susceptible to noise-related effects on calculated emissions mass at ultra-low concentrations.</p>	 <p>NOx Concentration (ppm)</p> <p>Time (s)</p> <p>Average concentration: ~0.5 ppm</p>
Raw Gas	A high-concentration peak occurs initially, followed by low concentrations after catalyst warm-up. A wide measurement range including high concentrations is required. Range switching is not permitted during regulatory testing.	
	 <p>NOx Concentration (ppm)</p> <p>Time (s)</p> <p>Average concentration: ~10 ppm</p>	 <p>NOx Concentration (ppm)</p> <p>Time (s)</p> <p>Average concentration: ~5 ppm</p>

High-Sensitivity CLD

CLD : Chemi-luminescence Detection

- Sample gas (NO) and O₃ gas is introduced into the reaction chamber at the same time.

NO is oxidized by O₃ to form NO₂.



10% of the NO₂ is in excited state.

NO₂^{*} return to ground state by emitting a photon.




- NO₂ to NO Converter


NO₂ is reductive to form NO to be detected by CLD.

This is not a catalyst, material consuming itself.


NOx converter efficiency check needs to be performed regularly.



Chemi-luminescence Wavelength



Structure of NO detector




To get more Chemi-luminescence:

- Optimal Sample Gas Flow
- Optimal O₃ Flow
- Reaction Chamber Temp.
- Reaction Chamber Press.
- Electrical Processing (Wide Range)

➔ High-Sensitivity CLD

Standard CLD vs. High-Sensitivity CLD



Zero Noise Performance



Compared with the standard NOx analyzer, the low-concentration NOx analyzer has a lower zero-point noise level.

Measurement Using High-Sensitivity CLD

Zero Stability and Impact on Mass Emission



Effect on Emission Mass

Standard NOx analyzer: 6%^{*2}

Low-concentration NOx analyzer: 0.3%^{*2}

^{*2} Impact relative to the 15 mg/mile regulatory limit


► The low-concentration NOx analyzer greatly improves the repeatability of emissions-mass measurement.

^{*1} LA#4 (UDDS): Los Angeles Urban Dynamometer Driving Schedule
A dynamometer driving mode simulating urban driving in Los Angeles

System Overview for 0.5mg/mile U.S. Tier 4 PM Standard

Combination of several technical optimization is possibly required

Single PM filter over multiple-phase tests



Tunnel temp. management with CVS




Increased face velocity (140cm/s) with advanced Dilution Sampler

Filter Face Velocity	Phase-1: 43%	Phase-2: 100%	Phase-3: 57%
Low FFV	43 cm/sec = 27 lpm	100 cm/sec = 63 lpm	57 cm/sec = 36 lpm
High FFV	60 cm/sec = 38 lpm	140 cm/sec = 88 lpm	80 cm/sec = 50 lpm
	Min FFV	Max FFV	

Tier 4 Upgrade for PM 0.5mg/mile



Tunnel flow rate flexibility with CVS



Automated filter weighing with Robot weighing system



Measurement of Low-Mass PM

Ensure sufficient PM collection and minimize variations in the measurement and weighing systems


Classification	Control Point
 Dilution Ratio Control	Maintain sufficient PM mass on the filter by using a low dilution ratio.
 PM Sampling Method	Use a single filter for three-phase sampling.
 Flow Rate and FFV Control	Combine flow weighting with a high filter face velocity (FFV).
 Background PM Control	Measure tunnel background PM before and after each test.
 Exhaust Line and Flow-Path Control	Minimize PM losses, condensation, and electrostatic effects in the exhaust line.
 Filter Selection	PTFE, reducing EC/OC effect but careful static-charge required. TX40 is less static-charge but EC/OC background effect must be considered.
 Weighing Accuracy Control	Reduce weighing errors by using the ABBA method and automated weighing.
 Environmental and Cleanliness Control	Continuously control temperature, humidity, dew point, and cleanliness.

CO₂ Measurement Challenges

Low CO₂ Concentration with High Moisture

Characteristics of Low-Carbon and Carbon-Neutral Fuels


Increased moisture content in exhaust gases



Comparison of Moisture Content in Exhaust Gases Between Fossil Fuels and Low-Carbon and Carbon-Neutral Fuels

Characteristics of PHEVs

Decreased CO₂ concentration in the CVS bag




UDDS : Urban Dynamometer Driving Schedule
(Driving cycle on a dynamometer simulating urban driving)


Comparison of exhaust emissions between a conventional vehicle and a PHEV

CO₂ Measurement Challenges (NDIR)

Water Interference



Conventional vehicle exhaust





Hybrid vehicle x low-carbon / decarbonized-fuel exhaust

- CO₂ concentration in the bag decreases while the relative proportion of water increases, which makes the measurement more susceptible to water interference.
- Therefore, in BAG measurements without sample dehumidification, low-concentration CO₂ is relatively susceptible to water interference.

* In direct measurements with sample-gas dehumidification, water-interference effects are not an issue.

Measurement Using a Moisture Compensated CO₂ Analyzer

Compensation Function



Structure of an NDIR (Non-dispersive Infrared) analyzer with compensation detector

Summary

New Measurement Technology and Response Strategies for Euro 7 and U.S. Tier 4

Item	Summary
Background	As emissions requirements become increasingly stringent, the following new measurement requirements are emerging:
1. Ultra-low NO _x measurement	Provide excellent zero-point stability, high sensitivity, low noise, and a wide dynamic range for ppb-level measurements.
2. Low-mass PM measurement	Ensure sufficient PM collection and minimize variations in the measurement and weighing systems to support measurements at the 0.5 mg/mile level.
3. Low-CO ₂ with high moisture measurement	Use an analyzer with water-interference compensation to address lower CO ₂ concentrations and increased water content in exhaust gas.

Omoshiro-okashiku
Joy and Fun

眞峰

Thank you

Cảm ơn

감사합니다

ありがとうございました

Dziękuję

धन्यवाद

Grazie

Merci

谢谢

ආචාර්ය

நன்றி

Gracias

Obrigado

Σας ευχαριστούμε

Děkuji

Teşekkürler

شکرا

Tack ska ni ha

Danke

Большое спасибо